

High Throughput Analysis of Food and Oil Contaminants


Fully automated
analysis of MCPD,
2-MCPD, 3-MCPD and
Glycidyl fatty esters.





Analysis of Contaminants designed for your Challenges


Process contaminants of the food industry were evaluated by the EFSA (European Food Safety Authority). Zinsser Analytic invented a new system for contaminant analysis in food and oil samples in collaboration with a food quality and safety lab, addressing your everyday workflow challenges with the following substances:

 **2-Monochloropropan-1,3-diol (2-MCPD)**


 **3-Monochloropropan-1,2-diol (3-MCPD)**

 **Esters of 2-MCPD and 3-MCPD (2-MCPD-E, 3-MCPD-E)**


 **Glycidyl esters (GE)**




2-MCPD




2-MCPD-E



3-MCPD



3-MCPD-E



GE

Those substances are common food contaminants in the following oils and processed foods:



Palm oil



Sunflower oil



Fine bakery wares



Corn oil



Soybean oil



Smoked fish and meat products



Grape seed oil



Avocado oil etc.



Fried or roast meat



Coconut fat



Margarine



Chips, crisps and dough based analogues



Olive oils



Bread and rolls



Infant and follow-on formulae

With our system you will have a lot of advantages combined:


✓ Reproducible and reliable results

✓ Certified workflow

✓ A larger throughput

✓ Less health impaired staff from manual pipetting


✓ Simplified sample handling




INTEGRATION

- 1 Left robotic arm with standard tips and pneumatic vial gripper.
- 2 Wash station for pipetting tips.
- 3 Reagents rack for 24 vials (40 mL or 80 mL).
- 4 DESYRE-Mix® racks for 24 vials each.
- 5 Lidding station for Desyre blocks.
- 6 Ultrasonic bath with 24 positions for vials.
- 7 DESYRE-Mix® for shaking of racks.
- 8 Workbench cut-out for centrifuge racks.
- 9 Centrifuge racks with six vial positions each.
- 10 Right robotic arm with gripper for lids and racks.
- 11 Stacker for lids.


CERTIFIED AUTOMATABLE PROCESSES



Ultrasonic bath



Centrifuge racks



Lidding station

The presented automated high throughput sample liquid handling system performs four different well established methods to analyze 3-MCPD and 2-MCPD (free and ester form), as well as glycidyl ester in oil and fat containing samples.

With the standard pipetting tips almost every liquid can be handled. The gripper for racks and vials can handle every needed vial and rack type, as well as rack lids. Sonication of vials could be performed with the incorporated ultrasonic bath. The needed centrifugation steps are performed with the centrifuge underneath the workbench with a corresponding cut-out for the gripper to install the racks with the centrifuge.

The methods applied can be changed, depending on the certified process desired: Alkaline transesterification of glycidol to 3-MCPD (DGF method, Cd 29c-13), alkaline transesterification of glycidol to 3-MBPD (SGS 3-in-1 method, Cd 29b-13), acidic transesterification (Unilever combined method, Cd 29a-13) and single compound analysis.

SCHEME 1

Unilever combined method (Cd 29a-13) for the simultaneous quantification analysis of 2- and 3-MCPD and GE in edible oils and fats.

A: PREPARATION

1. Sample aliquot dissolved in THF
2. Add internal standards
3. Add acidified solution of NaBr (0.1%, w/v)
4. Heat to 50 °C and vortex for 15 min
5. Transfer organic phase into new vial
6. Evaporate organic solvent



B: ACIDIC TRANSESTERIFICATION

1. Reconstitute in THF
2. Add H₂SO₄/methanol (1.8%, v/v)
3. Heat to 40 °C and vortex for 16 h
4. Neutralize solution
5. Add NaSO₄ solution (20%, w/v)
6. Discard aqueous phase



C: DERIVATIZATION IN ORGANIC PHASE

1. Add phenylboronic acid in diethylether
2. Incubate for 5 min at room temperature in ultrasonic bath
3. Evaporate organic phase
4. Reconstitution in iso-octane
5. GC-MS analysis

SCHEME 2

SGS 3-in-1 method (Cd 29b-13) for the quantification analysis of 2- and 3-MCPD esters and GE in edible oil and fats.

Assay A

A: ALKALINE TRANSESTERIFICATION

1. Dissolve sample in diethylether
2. Add internal standards
3. Mixture is cooled at -22 °C to -25 °C
4. Add methanol/NaOH (0.25%, w/v)
5. Vortex 16 h at -22 °C to -25 °C
6. Add acidified NaBr solution (60%, w/v)



B: EXTRACTION

1. Add iso-hexane
2. Vortex
3. Discard organic phase
4. Add diethylether/ethylacetate (3/2, v/v)
5. Vortex
6. Discard aqueous phase



C: DERIVATIZATION IN ORGANIC PHASE

1. Add phenylboronic acid in diethylether
2. Incubate for 5 min at room temperature in ultrasonic bath
3. Evaporate organic phase
4. Reconstitution in iso-octane
5. GC-MS analysis

Assay B

A: ALKALINE TRANSESTERIFICATION

1. Dissolve sample in diethylether
2. Add internal standards (different from assay A)
3. Mixture is cooled at -22 °C to -25 °C
4. Add methanol/NaOH (0.25%, w/v)
5. Vortex 16 h at -22 °C to -25 °C
6. Add acidified NaBr solution (60%, w/v)



B: EXTRACTION

1. Add iso-hexane
2. Vortex
3. Discard organic phase
4. Add diethylether/ethylacetate (3/2, v/v)
5. Vortex
6. Discard aqueous phase




C: DERIVATIZATION IN ORGANIC PHASE

1. Add phenylboronic acid in diethylether
2. Incubate for 5 min at room temperature in ultrasonic bath
3. Evaporate organic phase
4. Reconstitution in iso-octane
5. GC-MS analysis

SCHEME 3

DGF method (Cd 29c-13) for the differential quantification of 3-MCPD esters and GE in edible oils and fats.




SOFTWARE TACKLING EVERY CHALLENGE IN YOUR PROCESS

Every process step is designed modularly. Therefore a variety of processes, using differing chemicals and solvents, may be carried out. The customized sequential arrangement of the process is composed with a specified MS Excel® file or our Worklist-Generator which is executed with our WinLissy®-software.



Worklist-Generator

With the Worklist-Generator you can create your individual processes with many adjustable parameters for sub steps. Drag and drop the needed methods to the desired positions of the worklist.




Method Definition

Input the desired time or anything the method requires, while you are defining the process.



Pipetting-Editor

For pipetting steps, the Pipetting-Editor is opened for the definition of racks and number of steps. Every filled position is indicated blue and the selected positions are highlighted green.



WinLissy®-Software

With our WinLissy®-software every prepared worklist of the Worklist-Generator is executable. Communication with customer specific databases or LIMS is easily conceivable.

GENERAL CONTACT

EMEA

Zinsser Analytic GmbH
Schwalbacher Str. 62
65760 Eschborn
Germany

Tel: +49 6196 58693 0
Fax: +49 6196 58693 80
info.zinsser-analytic@gardnerdenver.com

AMERICAS

Zinsser NA, Inc.
19145 Parthenia Street
Suite C
Northridge, CA 91324
USA

Tel: +1 818 341 2906
Fax: +1 818 341 2927
info@zinsserna.com

HEADQUARTERS

Gardner Denver Thomas GmbH
Livry-Gargan-Str. 10
82256 Fürstenfeldbruck
Germany

Tel: +49 8141 2280 0
Fax: +49 8141 8892136
thomas.de@gardnerdenver.com

Please check out all our brands for your mission-critical flow control technologies:



OEM Pumps
www.gd-thomas.com



Vacuum Pumps & Systems
www.welchvacuum.com



Syringe Pumps
www.tricontinent.com



Labautomation
www.zinsser-analytic.com



Syringes
www.microsyringes.com